



GUIDELINES FOR LIVE DEMONSTRATIONS (LIVE DEMO) ON HUMAN MODELS

Live Demonstrations: RA, PoCUS, Chronic Pain

Definition:

These sessions are **NOT HANDS-ON** Workshops!

Delegates will interact with the Demonstrator(s) but **NOT** scan themselves! These Demonstrations will take place in a **LARGE** room with **GOOD** visibility to the screen (projected pictures from the Ultrasound machine and the hand of the Demonstrator while moving) and to the Demonstrator / Human Model.

This session is the **Combination of a Lecture with Live Scanning** by the Demonstrator(s) and Interaction with the Audience. Therefore, they are very attractive for delegates, leading usually to overcrowded sessions.

Demonstrators: One or two Demonstrator(s) per Live Demonstration on Human Models Session.

Chairman: There will be **NO Chairman** for these sessions. The Demonstrators will guide the session, scan, explain and keep the time themselves. In the case of two demonstrators, they will share the work according to their preparation, and as such, in case of 2, it is advised they get in touch with each other.

Material needed:

LARGE lecture room! This has already been taken care of by the organizers.

- Projection from Ultrasound machine (US) to screen (showing the hand moving as a window in the projection screen) and from laptop to screen.

In every Live Demonstration there has to be:

- **One** male or female model with **one HIGH END US** (meaning the best possible US from the corresponding company with high resolution as the pictures have to be projected to a screen).
- Complete US equipment (battery loaded, working connection for power supply), enough gel for all demonstrations, appropriate towels to clean the US-probe.
- Stretcher or other comfortable surface and pillows for positioning the models for a duration of up to 50 minutes of scanning in different positions per live demonstrations.
- Enough towels and blankets to cover the models.
- At least one pillow per stretcher or massage bed / model.
- One bathing gown per model (for the time between the live demonstrations, with the models being told to wear easily removable sport wear).
- Enough paper towels to clean the gel from the models for all live demonstrations.

Per Live Demonstration:

- Water-soluble skin pen to paint on the models, different needles and catheters for demonstration (supplied by the companies at the congress)
- 2/3 chairs per live demonstration (for instructor and for model between the live demonstrations)

**6TH WORLD CONGRESS
ON REGIONAL ANAESTHESIA
& PAIN MEDICINE**

Joint meeting with
**THE 40TH ESRA
ANNUAL CONGRESS**
esraworld2023.com



Audio – Visual support for all live demonstrations:

- There **MUST BE** projection from the US machine and the demonstrator's hand (as a window on the screen) and from a laptop (need for beamer, laptop on place to upload a presentation or to connect a USB). Also, a connection to use the own laptop or Mac of the demonstrator is needed for a short presentations during the scanning (there will be changes from presentations from the laptop to the US and vice versa).
- Two headset microphones needed per live demonstration
- Two laser pointers (ideally green and not red light) to show details on the screen
- Microphones and hostesses for interaction with the audience

Perfect set-up: See document «Live demonstration set-up»

Support from US companies: ONLY HIGH END US for best possible imagine resolution during projection

For the **PAIN LIVE DEMONSTRATIONS** consider the following companies: GE, SAMSUNG, TOSHIBA, SONOSITE

At least 1 industry person per Live Demonstration present per US-label/brand used (Sonosite, BK, GE, Mindray etc.) until the instructor feels comfortable with the US

Projection for Presentations: NEEDED. See above

Computer: NEEDED. A Computer / laptop to install a presentation or a USB and the possibility for connection to an own laptop (Mac). See above.

Responsible Persons on Site

- **Program Coordinator: Alain Delbos**, Chair of Scientific Committee

Cell Phone Number: +33 670795088

- **Scientific Committee: Jose Aguirre**, Chair of Hands – On Clinical WS

Cell Phone Number: +41 76 331 11 48)

- **PCO: MCO Congres Group – Viviane Barbarisi** - Cell Phone Number: +33 670069826

• **Back - Up:** A back – up to possibly cover unplanned absences is more difficult than for Clinical Workshops, due to the special expertise needed and the duration of the Live Demonstration. Therefore, some of the live demonstrations are scheduled with two demonstrators.

However, the PCO is advised to be in touch with the faculty at the beginning of the congress to check if the demonstrator is available.

Who to inform in the case of problems?

To avoid delays and other problems, ALL information concerning:

- Problems with models / ultrasounds
- Non-show by demonstrators
- Cancellation by demonstrators prior to the workshop
- Problem with rooms which cannot be solved prior to a workshop
- Problems with the number of candidates (Workshop is over-booked, need for more instructors etc) Has to be shared with the following on behalf of ESRA

Jose Aguirre: +41 76 331 11 48

Kalliopi Platsi: +30 6945026525

PCO (MCO Congres Group) – Viviane Barbarisi - Cell Phone Number: +33 670069826